

# FASCIAL TRACTION

## Fascial traction system | Mod. 8540

- Symmetrical and controlled abdominal closure.
- Allows traction at any angle.
- Regulated (marked) bars and traction system for traction control.
- Plate with different holes for suture support.
- Full plate movement thanks to the ball joint.
- Precise traction: rack-and-gear mechanism.
- Modular system for storage.



	UDS.	REF.	DESCRIPTION
<b>A</b>	2	B230	Radial inclined flange
<b>B</b>	2	AR160	Bar adapter d16x300 mm
<b>C</b>	2	R210	Bar d16x550 mm
<b>D</b>	2	P840	Fascial traction system
<b>E</b>	2	B416	Toothed flange d16-d16
<b>F</b>	1	R250D	Detachable bar d16x750 mm

PRESENTATION  
VIDEO



# FASCIAL TRACTION

## LAPAROSCOPIC SURGERY KIT

DESCRIPTION	QTY	CODE
X-ONE Hasson Trocar 12mm Standard Cannula with Balloon	1 Trocars	300-HB12
X-ONE Optical Trocar 11mm Standard Cannula	1 Trocars	300-V11
X-ONE Optical Trocar 5mm Standard Cannula	2 Trocars	300-V5
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 37mm (HR) - USP 1 EP 4 - 30cm</i>	3 Threads	FQ467HHAE
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 26mm (HR) - USP 2/0 EP 3 - 20cm</i>	3 Threads	FS412HHAC
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 37mm (HR) - USP 0 EP 3,5 - 45cm</i>	3 Threads	FR465HHAE
ASSUFIL - Absorbable, braided polymer of glycolid acid (60-90 days) <i>Taper J Needle 30mm (JR) - USP 1 EP 4 - 60cm</i>	3 Threads	FR82A AM
ASSUSILK - Non-Absorbable Braided Silk Suture. <i>Straight Cutting Edge 60mm (GS) - USP 1 EP 4 - 3x100cm</i>	4 Packs	03KR719 AV
ASSUMESH - LOW DENSITY, Oval Non-Absorbable monofilament polypropylene 55gr/m <sup>2</sup> Mesh - 30x40cm	1 Mesh	AM30420

## OPEN SURGERY KIT

DESCRIPTION	QTY	CODE
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 37mm (HR) - USP 1 EP 4 - 30cm</i>	3 Threads	FQ467HHAE
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 26mm (HR) - USP 2/0 EP 3 - 20cm</i>	3 Threads	FS412HHAC
FILBLOC - Absorbable, Barbed Suture in Polidioxanone unidirectional with final lock system (180-210 days) <i>Taper 1/2 C. 37mm (HR) - USP 0 EP 3,5 - 45cm</i>	3 Threads	FR465HHAE
ASSUSILK - Non-Absorbable Braided Silk Suture. <i>Heavy Taper 1/2 C. 48mm (HRG) - USP 1 EP 4 - 3x100cm</i>	4 Packs	03KR479 AV
ASSUMESH - LOW DENSITY, Oval Non-Absorbable monofilament polypropylene 55gr/m <sup>2</sup> Mesh - 30x40cm	1 Mesh	AM30420